# COMPOSITE AUTO

## CARACTERÍSTICAS PRINCIPALES:

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# DENTALFLUX

Scientific Pharmaceuticals  
3221 Producer Way Pomona  
CA USA

Quality Firts International Ltd  
20 Eversley Rd Bexhill-On-Sea  
East Sussex Inglaterra TN401HE UK

## OUTSTANDING FEATURES OF THE MATERIAL

High filler content	Low pulp irritation potential
High mechanical strength	Excellent X Ray opacity
High wear resistance	Good polishability
Long Selg life	Resistance to staining
Excellent marginal adaptation	Excellent aesthetics
Low contraction	Easy handling

## PROPERTIES OF CURED RESTORATIVE

	Requirements of ADA N°27	Sci Pharm composite
Flexuran strength	50 Mpa	120 Mpa
Hardness (Barcol)	Not specified	> 98
Coefficient of thermal expansion	Not specified	$31 \times 10^{-6} \text{ cm}^2/\text{C}^\circ$
Opacity/translucency Factor C70	0,35-0,75	0,45
Water absorption	max 0,7 µg/cm <sup>2</sup>	0,5 mg/cm <sup>2</sup>
X-Ray opacity	> twice the X ray opacity of Al	Complies
Color stability	Discoloration perceptible with Difficulty	Complies
Working time	1,5 min	2 min
Setting time	Max 8 min	4-5 min

## GENERAL INFORMATION

Paste paste composite self cure has the following main characteristics: high proportion of filler/resin, low water absorption, high mechanical strength and hardness.

### ADDITIONAL TECHNICAL INFORMATION REQUIRED BY THE ADA REVISED SPECIFICATION

- 1.- Principle organic component of material Bis-GMA and polyglycol dimethacrylate resin blend
- 2.- Particle sizes of inorganic filler: 100% below 30 microns, 90% below 10 microns, 50% below 2 microns
- 3.- Volume % of filler in the restorative material 63%

### CLINICAL PROCEDURE

Cavities are prepared in the conventional manner. In deep cavities it is recommended to place a cavity liner calcium hydroxide based or Fluoroseal. Etch the margins for an optimal adherence. Before placing the composite an adhesive must be placed. If a+b adhesive is used, a small quantity of parts A and part B adhesive must be mixed and placed in the cavity

### SPECIAL RECOMMENDATIONS

**CLASS 1:** Slightly overfill the cavity. Place a plastic strip over the restoration and instruct the patient to bite down for one minute. Remove excess restorative and trim the flash. The restoration will be ready to finish four minutes after inserting the paste into the cavity.

**CLASS 2:** the use of plastic strip is recommended for interproximal separation and as a matrix. Compress the restorative in the cavity with a strip. Hold the strip three minutes until the restorative sets. Trim and finish after four minutes.

**CLASS IV:** The use of bonding agent over etched enamel is a necessary step in this kind of restoration. To provide adequate bonding strength, a relative large area of enamel (2-3mm over) should be etched and primed with bonding agent. Retention may also be increased making small undercut on the lingual side of the tooth. In other situations, the use of pins placed in the dentin may be necessary even crowns.

**MIXING AND APPLICATION INSTRUCTIONS:** With the opposite ends of a spatula, dispense equal quantities of part A and B, equal parts A and B results in universal color. The proportion of part A and B can be modified slightly to adjust the color. Spatulate during 20 to 30 seconds obtain a mix of uniform color. Place the composite. The working time is 2 min at 23°C. If the material is cooler, the working time is longer. If matrix bands are used, they should be not removed before four minutes. Crowns must be not removed before six minutes. The polish usually is done after five minutes.

**STORING:** keep the storage below 23°C. Refrigerate when the material is not in use. Be careful not to leave outside refrigeration 1 summer time.

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