EUGENOL

(ES)

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(EN)

INDICATIONS:

When eugenol is mixed with pure zinc oxide or any zinc oxide based cement, a dental cement is obtained.

These cements are used as endodontic fillers and other various uses depending of the formulation of the powder, that better water resistance and hardness.

INSTRUCTIONS FOR USE:

Usually the proportions of the mixture is 3:1- 6:1 powder eugenol, for ZOE cements, but with the high variety of formulations marketed, the recommendations of the powder manufacturer must be followed.

ACTION:

The eugenol reacts with the zinc oxide following this reaction:

- a) $ZnO + H_2O \rightarrow Zn(OH)_2$
- b) $Zn(OH)_2 + 2 \text{ Eugenol} \rightarrow Zn \text{ Eugenate} + H_2O$

This is a common acid/alkaline reaction but not so clear. It does not liberate heat, it is so slow, needs several hours to complete. The result is a non resistant matrix where the salt of eugenate acts a matrix and the remaining zinc oxide as filler. The remaining eugenol causes inhibition in polymerization processes of resins. And as can be observed water is essential to start the reaction.

PRECAUTIONS:

The rests of zinc oxide eugenol cements inhibits resins polymerization.

COMPOSITION: EUGENOL 99% min.

Storage: Keep the bottle closed and under 25°C.

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