FLUROSEAL

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(EN)

a) Chemical composition:

Fluroseal does not contains calcium hydroxide $Ca(OH)_2$, but a calcium salt of apatite hydroxide ($Ca_{10}(PO_4)_6(OH)_2$). This compound has a biocompatibility long time tested in implantology procedures. The matrix of the formulation is Diuretane dimetacrilate, which provides the cured product with a hardness higher than the glass ionomers, as well with an excellent chemical union with composites. Another advantage of this matrix is the great thermal isolation. This thermal isolation is very important to avoid the heat produced in resins (composites) curing, reaches the vital pulp. Also a great resistance against chemical products is offered. It can be seen when it is exposed to etchers and no effect is observed. Fluroseal contains sodium fluoride. With the time frees the fluor inside the teeth. This liberation reduces the sensibility due to a harder effect in the dentin. Fluroseal is radiopaque for future RX texting.

b) Pulp protection:

The dentin generation is produced when a strongly alkaline irritation, and Ca^{++} ions are available. The Apatite hydroxide calcium salt ($Ca_{10}(PO_4)_6(OH)_2$), is chemically closed with the well known calcium hydroxide ($Ca(OH)_2$), also produces a strong alkaline medium, and frees Calcium ions. Anyway the secondary dentin formation has been observed clinically.

Hardness	Compresion (PSI)	Tensil Strenght (PSI)	Solubility
FLUROSEAL	20.000-20.500	4.000-5.000	0.05%

INSTRUCTIONS:

1. Prepare the cavity as usual. Is important a dry cavity to assure an optimal adherence.
2. If the pulp is so close, a little quantity of calcium hydroxide should be placed under the fluroseal. Anyway place a thin layer (0,5 mm) of **FLUROSEAL** and cure for 30 seconds.
3. Place a new layer of fluroseal thicker, and cure it.
4. Continue with the usual way.

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